



# CLAIND LABORATORY GAS GENERATORS

## nitrogen for LC-MS

The important considerations of safety, convenience and cost savings influenced the design of the new nitrogen generator tailored for LC-MS analysers.

The generator is not hazardous and can be installed in the laboratory if required, although it is normally sited remotely with a nylon pipe feed (6 mm dia.) to the instrument.

The advanced technology developed by CLAIND produces pure nitrogen, allowing the **N<sub>2</sub> LCMS** generator to supply **any LC-MS** produced by **any manufacturer!**

When compared to costs of bought-in nitrogen supplies the **N<sub>2</sub> LCMS** should show a capital repayment time **less than 16 months**. Thereafter, the nitrogen is supplied virtually free of charge.

**Independence** is the main feature of these generators. They only require a supply of mains electricity for automatic, unattended operation. Pressing a button is the only operation to start up the generator, which operates automatically 24 hours a day, 365 days a year.

Maintenance is limited to changing a few filter elements taking less than 30 minutes per year. The system regenerates itself automatically giving you a continuous, uninterrupted supply of gas.

### MAIN FEATURES

- Ready to use just few minutes after switching on
- Gas production proportional to the consumption with automatic 'Stand-by' when demand falls
- Compact instrument, just 40 cm wide supplied **with or without oil free air compressor fitted internally**
- Very low noise level
- 50 liters internal air buffer included with the generator
- Direct connection to the users LC-MS



### TESTED WITH

WATERS

APPLERA/APPLIED BIOSYSTEM

VARIAN

THERMOFINNIGAN

MICROMASS

AGILENT

SHIMADZU

### WORKING PRINCIPLE

Thanks to the **PSA (Pressure Swing Adsorption)** principle the generator produces nitrogen by compressing ambient air and passing it into a Carbon Molecular Sieve bed (CMS).

Inside the CMS bed, oxygen, moisture, CO<sub>2</sub> and other 'contaminants' are trapped allowing nitrogen to pass through into a holding reservoir.

Nitrogen from the holding reservoir is regulated to a fixed flow and pressure before exiting the generator.

**Claind is one of the worldwide major manufacturer of the PSA technology.**

# LC-MS *series*

## gas characteristics

Nitrogen purity according to the LC-MS manufacturers specifications\*

Outlet pressure: 7 barg

\* LC-MS generators have been tested from Claind with all the main LC-MS now on the market

## Models available

**N<sub>2</sub> LC-MS I**: nitrogen generator with built in compressor

**N<sub>2</sub> LC-MS 0\***: nitrogen generator with external air supply

\* suggested in case of 2xLCMS or in case of APCI interface when higher flow is requested

## technical specifications

Technical specifications	N <sub>2</sub> LC-MS I	N <sub>2</sub> LC-MS 0
<b>FLOW RATE</b>	max 15 NI/min	max 38 NI/min
<b>POWER RATING</b>	850 VA	80 VA
<b>ELECTRICAL SUPPLY</b>	230 Vac - 50 Hz *	230 Vac - 50 Hz *

\* 115 Vac - 60 Hz available on demand

Noise level: <60 dB

Operating temperature: between 5°C and 40°C

Safety protection: > 11,5 bar

Consumables: filters and silencers

### AIR SUPPLY CHARACTERISTICS\*

Inlet Air flow: .....min 110 NI/min

Inlet Air pressure: .....min 8,5 / max 10 barg

Dew point: .....< +3° C

Particles: .....< 0,01 µm

Oil vapors: .....< 0,01 mg/m<sup>3</sup>

\* Only for N2 LC-MS 0

## dimensions

	N <sub>2</sub> LC-MS I	N <sub>2</sub> LC-MS 0
<b>Height</b>	119 cm	119 cm
<b>Width</b>	40 cm	40 cm
<b>Depth</b>	80 cm	80 cm
<b>Weight</b>	119 kg	112 kg

### CLAIND s.r.l.

Via Regina 24, 22016 Lenno - Italy Tel: +39 (0) 34456603 Fax: +39 (0) 34456627 e-mail: [info@claind.it](mailto:info@claind.it) - [www.claind.it](http://www.claind.it)

Claind is a trademark of Claind s.r.l.. Copyright: Claind limited 2003. Document version: PSLCMTEN200307-2. The company reserves the right to change the above specifications according to the product development. For updated literature please refer to Claind sales representatives.